

G350 Fully Automatic Filter Bag Fat Extractor

G350 — Fully Automatic Filter Bag Fat Extractor (Soxhlet Alternative)

Overview

Fully automatic filter-bag fat extractor with automated solvent adding — covering the complete process of sealing, soaking, extraction and solvent recovery, for up to 15 samples per batch.

Features

- Automated solvent adding.
- Full process: sealing, soaking, extraction and solvent recovery.
- User-defined schemes: 250 sets.
- Repeatability RSD ≤1.0%.

Specifications

Measuring range	0.1-100%
Capacity per batch	Up to 15
Sample size	0.1-5 g
Temperature range	Ambient-105 °C
Solvent recovery	≥90%
Solvent adding	Automated
User-defined schemes	250 sets
Repeatability	RSD ≤1.0%

All technical parameters are subject to change without notice.

Application fields

Total fat / crude fat determination in food, feed, cereals, meat and oilseeds by automated solvent extraction.

Specifications as published at www.sino-sonnen.com. Final configuration is confirmed per order — request a quotation: info@sino-sonnen.com · R&D: Jiading District, Shanghai · Manufacturing: Luoyang, Henan, China