

Gr08 Automated Graphite Block Digester

Gr08 Automated Graphite Block Digester

Overview

This kind of graphite digester is adopted with the economic high purity graphite heating method and the anti-oxidation processing. Comparing to that of the aluminum block, the price of the graphite digester is more acceptable to most of the laboratories.

Features

- Heat insulation methods: Ceramic insulation and air duct insulation.
- Heating method: Infrared radiation and graphite conduction heating.
- Temperature is adjustable and constant.
- Alarm of over voltage, over current, over heating, over temperature, and auto malfunction alarm.
- Accurate PID outlet temperature control and safe over temperature protection are specially designed.
- The outer surface is painted with teflon coating to prevent the apparatus from being corroded by the acid gas or liquid generated during the digestion.
- The exhaust manifolds should be connected to the water aspirator supplied, or to a suitable scrubber.
- Sonnen's E Series Off-gas Absorption Apparatus: connection of the digester and off-gas absorption system, to neutralize and absorb the acid fume generated during digestion.

Specifications

Temperature range	Room temperature to 450°C
Temperature stability	±1°C
Heating method	High purity graphite conduction heating
Thermal insulation material	Ceramic fiber
Digestion tube volume	300 mL
Processing capacity	8 samples per batch
Power source	AC220VAC±10% 50Hz
Rated power	2000W
Dimensions(mm)	340×320×165
Net weight	18 kg

All technical parameters are subject to change without notice.

Application fields

This graphite digester is widely applied in the industries of food, pharmaceutical, agriculture, forestry, environment protection, chemical engineering, biochemistry and etc. For those colleges and universities, and other scientific research departments, which will handle the sample digestion before chemical analysis of the soil, fodder, plant, seed, mine and etc, it is better to equip Sonnen's Kjeldahl apparatus together with this graphite digester.

Specifications as published at www.sino-sonnen.com. Final configuration is confirmed per order — request a quotation: info@sino-sonnen.com · R&D: Jiading District, Shanghai · Manufacturing: Luoyang, Henan, China