

T11 Intelligent Multifunctional Titrator

T11 Intelligent Multifunctional Titrator

Overview

Sonnen's titrator focus on precision, intelligent and modularization, Titration unit achieve "build blocks", In functionality different card for client's choice, contain potentiometric, PH, color and other functions.

Features

- Building Block Design, T01 titron, T11E Electric Potential board, T11C color card, T02 magnetic stirrer can be willfully combinate or individually work. Concerning the separate use, please refer to the details of T03(which contains T01 & T02).
- High brightness color touch screen(480 X 320 pixels), individual operation interface, can set 12 graphical shortcut keys.
- English and Chinese operation interface are alternative, Print interface reserved.
- Intuitional"Building Block"flexible combination, operator can select the relative titration station and function according to actual need, and upgrade anytime.
- Built-in multiple analysis method.
- Host bits can control twelve titron of T03 by bluetooth.
- Host can connect to PC by WLAN.
- The titration unit of T03 can work individually, mouse control titration.
- The alternative communication interfaces: two USB interfaces, two RS232 interfaces, CAN bus interface, WLAN interface.
- Different combination meets different clients' demand.

Specifications

The full specification sheet for this model is confirmed per order — contact us for the detailed configuration.

Application fields

Chemical analysis laboratory of various industries

Specifications as published at www.sino-sonnen.com. Final configuration is confirmed per order — request a quotation: info@sino-sonnen.com · R&D: Jiading District, Shanghai · Manufacturing: Luoyang, Henan, China